

Measurement of the Radioactive Nuclides for the Medical Linear Accelerators based on the Spectroscopy using In-situ HPGe detector

Na Hye Kwon, Dong Wook Kim, Jinsung Kim, Young Jae Jang, Kum Bae Kim, Sang Hyoun Choi*


Department of Radiation Oncology, Yonsei University College of Medicine, Seoul, South Korea

Severance


Activation in Linear accelerator (LINAC)

- **Purpose**

- In the case of using X-ray energy above 8 MV, photoneutrons are generated by the photonuclear reaction
- The components of the linear accelerator (LINAC) are “**ACTIVATED**”
- We measured the radioactive nuclides using an in-situ High Purity Germanium (HPGe) detector for the real-disposed LINACs and components.









Activation



- **LINAC > 8 MV energy, Photo-nuclear reaction**
- **Neutrons react with “Target, Filter, Collimator, Jaw...”**
- **Components be activated**

Spectroscopy

- ✓ Spectroscopy (In-situ HPGe): The HPGe detector was positioned in front of the component and measured for 30 min for each component
- ✓ Targets for measurement: Disposed of LINACs (LINAC by Manufacturer: 1 Siemens, 2 Varian, 1 Elekta)

Manufacturer		Varian 1 (15 MV)			Siemens (15 MV)		
Component							
Dose rate [$\mu\text{Sv/hr}$]		15.6 \rightarrow 2.21 *5주 뒤)	4.10 \rightarrow 2.37	2.43 \rightarrow 2.49	26.9 \rightarrow 7.05	6.97 \rightarrow 3.64	7.65 \rightarrow 7
* Half life > 100 days		Target	15 MV Filter	Primary Collimator	Target	15 MV Filter	Primary Collimator
⁵¹ Cr	27.70 d						
⁵⁴ Mn	312.2 d						
⁵⁷ Co	270.9 d						
⁵⁷ Ni	35.65 h						
⁵⁸ Co	70.80 d						
⁵⁹ Fe	44.63 d						
⁶⁰ Co	5.271 y						
⁶⁴ Cu	12.70 h						
⁶⁵ Zn	244.26 d						
⁹⁹ Mo	66.02 h						
¹⁸¹ W	121.2 d						
¹⁸⁴ Re	38.00 d						
¹⁸⁷ W	23.83 h						
¹⁹⁶ Au	6.183 d						
¹⁹⁸ Au	64.80 d						

Spectroscopy

- ❖ The spectroscopy results for the LINAC “Varian 2” and “Elekta”, we couldn't figure the radio nuclides for the major components.
- ❖ LINAC Varian 2
 - Institution where disposed Varian 2, used the 6 MV for 2 months before the disposal.
 - Nuclides such as ^{60}Co or ^{54}Mn (long-lived radionuclide) was found from the lead block.
- ❖ LINAC Elekta
 - We measured only for 10 min for the each components. (lack of the cps)



Trans-SPEC-DX-100T



GC2018

Surface dose rate measurement



Thermo FH 40 G-L10
survey meter

- ✓ We measured the surface dose rate for the each components using a survey meter
 - Measured the dose rate at the surface
 - Measured at the day of disposal & 5 weeks after disposal

* Surface dose rate (Major component such as Target, Primary collimator, flattening filter for high energy)

Surface dose rate [μSv/hr]		Target	Primary Collimator		Flattening filter (for High energy)	
point of the measurement	Day of disposal	5 weeks after disposal	Day of disposal	5 weeks after disposal	Day of disposal	5 weeks after disposal
Elekta	3.318	0.459 ~ 0.578	27.1	0.152	0.161	0.128
Siemens	26.856	7.0488	7.65	7.002	6.966	3.636
Varian 1	15.5	2.38	4.59	3.46	4.1	2.69
Varian 2	0.209	-	0.153	-	127.764	-

* Background dose rate was around $\sim 0.125 \sim 0.135 \mu\text{Sv/hr}$

* For the other components, such as Ion chamber, MLC, Jaw, dose rate was around background level.

4. Conclusion

- **In this study, we measured immediate radionuclide analysis and dose rate evaluation using an In-situ HPGe detector and a survey meter for the medical LINACs.**
- **The difference in the level of radiation and the type of nuclide was confirmed depending on the period of use, energy, and workload.**
- **Therefore, in the process of relocating for the dismantling of the LINAC, it is considered to be safe exposure management when carried out 3 to 5 hours after stopping the operation.**
- **In future work, we will compare the nuclides acquired w/ Monte Carlo simulation.**

Acknowledgement: This study supported by the Nuclear Safety and Security Commission (No. 2003013) and by National Research Foundation of Korea (NRF) grant funded by the Korea government (Ministry of Sciences and ICT) (No. 2020M2D9A309417013).

Thanks for the attention !